

DLAT Monoclonal antibody

【货号】 AM20143

【规格】 50ul/100ul

【产品信息】

◆产品名称	DLAT Monoclonal antibody
◆规格	50ul/100ul
◆别名	1F2B9, Dihydrolipoyllysine-residue acetyltransferase component of pyruvate dehydrogenase complex, mitochondrial, DLTA, EC:2.3.1.12, ODP2
◆产品类型	Monoclonal
◆基因全称	dihydrolipoamide S-acetyltransferase
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human DLAT
◆理论分子量	70 kDa
◆实际分子量	70 kDa
◆Gene ID (NCBI)	1737
◆GenBank Accession Number	BC039084
◆UNIPROT	P10515
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:250-1:1000 for Immunohistochemistry (IHC) 1:200-1:800 for Immunofluorescence (IF/ICC)

【注意事项】 ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
◆本产品仅限科研使用。