

ALDH1A1 Monoclonal antibody

【货号】 AM20171

【规格】 50ul/100ul

【产品信息】

◆产品名称	ALDH1A1 Monoclonal antibody
◆规格	50ul/100ul
◆别名	1A10A2, 3-deoxyglucosone dehydrogenase, ALDC, Aldehyde dehydrogenase 1A1
◆产品类型	Monoclonal
◆基因全称	aldehyde dehydrogenase 1 family, member A1
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human ALDH1A1
◆理论分子量	55 kDa
◆实际分子量	55 kDa
◆Gene ID (NCBI)	216
◆GenBank Accession Number	BC001505
◆UNIPROT	P00352
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:25000 for western blotting (WB) 1:200-1:800 for Immunofluorescence (IF/ICC) 1:25-1:250 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。