

## PKM2-specific Monoclonal antibody

【货号】 AM20187

【规格】 50ul/100ul

### 【产品信息】

◆产品名称	PKM2-specific Monoclonal antibody
◆规格	50ul/100ul
◆别名	PKM2, PK-M2, PKM2 pyruvate kinase, muscle, PKM, Pyruvate kinase PKM
◆产品类型	Monoclonal
◆基因全称	PKM pyruvate kinase, muscle
◆经测试应用	WB, IHC, IF/ICC
◆宿主	Mouse
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human PKM
◆理论分子量	58 kDa
◆实际分子量	58 kDa
◆Gene ID (NCBI)	5315
◆GenBank Accession Number	NM_002654
◆UNIPROT	P14618
◆形式	Liquid
◆保存条件	Store at -20°C. Stable for one year after shipment.
◆建议稀释度	1:2500-1:15000 for western blotting (WB) 1:50-1:500 for Immunohistochemistry (IHC) 1:100-1:400 for Immunofluorescence (IF/ICC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
  - ◆本产品仅限科研使用。