

TMEM173/STING Polyclonal antibody

【货号】 AP20122

【规格】 50ul/100ul

【产品信息】

◆产品名称	TMEM173/STING Polyclonal antibody
◆规格	50ul/100ul
◆别名	STING, Endoplasmic reticulum interferon stimulator, ERIS, MITA, MPYS
◆产品类型	Polyclonal
◆基因全称	transmembrane protein 173
◆经测试应用	WB, IHC, IF/ICC, IF-P
◆宿主	Rabbit
◆种属反应性	Human, Mouse, Rat
◆免疫原	A synthesized peptide derived from human TMEM173
◆理论分子量	42 kDa
◆实际分子量	35-40 kDa
◆Gene ID (NCBI)	340061
◆GenBank Accession Number	BC047779
◆UNIPROT	Q86WV6
◆形式	Liquid
◆保存条件	Store at -20℃. Stable for one year after shipment.
◆建议稀释度	1:1000-1:5000 for western blotting (WB) 1:25-1:250 for Immunofluorescence (IF) 1:50-1:500 for Immunohistochemistry (IHC)

- 【注意事项】
- ◆为了您的安全和健康，请穿实验服并戴一次性手套操作；
 - ◆本产品仅限科研使用。